

DIAGNOSTIC PENLIGHT LED - LUXAMED

LUXAMED®

Ref. D1.211.112



 25/04/2025 21:43:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUXAMED®

Ref. D1.211.112

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

Features

MEASUREMENTS AND WEIGHT

DIMENSIONS AND WEIGHT			
*Subject to change without notice. For information purposes only, can be subject to change.			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 +33(0)4 42 90 31 31

 info@gsh-med.fr

 www.gsh-med.fr