

DIAGNOSTIC PENLIGHT LED - LUXAMED

LUXAMED®

Ref. D1.211.312



25/04/2025 21:36:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

DIAGNOSTIC PENLIGHT LED - LUXAMED

LUXAMED®

Ref. D1.211.312

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961697	HS Code	90189084
Warranty (months)	24	Country	Germany

MEASUREMENTS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
--------------------------	----	---------------------------	----

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

Colors Gray