

X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W

Ref. OTA1002



25/04/2025 21:37:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W



Ref. OTA1002

Description

These negatoscopes offer a mounting system for radios allowing simple installation and removal, and without consequence for radios. They can be placed on a piece of furniture or fixed on a wall. - The structure of the screen is made of high strength aluminum

- alloy 10 cm thick- The angles are protected by a steel plate, covered with ABS
- Methacrylic polymer screen- The body of the negatoscope is made of high strength aluminum alloy
- The inside and outside are covered with an antistatic sprayed powder
- Voltage / Frequency: 220V / 50-60Hz
- Luminance: 2,000 cd / m2
- Average lighting> 6300 lux
- Lighting life: 50,000 h
- Color temperature: 10,000 k- External dimensions: 80 x 56 x 10 cm
- Screen dimensions: 72 x 43.5 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040195865	HS Code	90105010
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	68	Bare Unit / Width (cm)	56
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	87	Bare Unit / Length (cm)	80
Packed unit / Weight	12	Bare Unit / Weight (kg)	8.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/04/2025 21:37:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr