

F.O. MACINTOSH BLADE WITH INTEGRATED LIGHT GUIDE, LUXAMED ®,

Ref. E1.414.012

LUXAMED®



28/04/2025 22:20:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LUXAMED®

Ref. E1.414.012

LUXAMED MACINTOSH laryngoscope blades with integrated fibre optics.

Description

Laryngoscope blade with integrated optical glass

- High-quality stainless steel construction for a long service life
- No dust collecting cracks – easy cleaning and sterilisation
- Maintenance-free
- 8.000 fibres to achieve the maximum service life and illumination
- Optimal light emission due to the largest possible cross section (3 x 7 mm)
- Brightness approx. 14.000 Lux at 3.7 V LED and 2.5 V LED
- Compatible with all handles according to ISO 7376 (green standard)
- Autoclavable up to 134°C – approx. 4.000 times
- Data-Matrix with most important information

Differents size available :

0 : E1.410.012

1 : E1.411.012

2: E1.412.012

3 : E1.413.012

4 : E1.414.012

5 : E1.415.012

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960546	HS Code	90189084
Warranty (months)	24		

MEASUREMENTS AND WEIGHT

For information purposes only. can be subject to change



GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr