

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

mediPlogs

Ref. WSAM_38BC



28/04/2025 22:54:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_38BC**

=====


Description


- Clogs suitable for modern hospitals, operating theatres and intensive care units.
- Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
 - The back of the heel is sloped and rounded to make walking easier
 - Ventilation system allows air to circulate
 - Toe prints to enhance hold and spurs to aid the circulation
 - 100% CFC-free polyurethane
 - Washable by hand or machine at 60°C
 - Sizes: 35 - 47

Features

| | | | |
|-----------------------------|---------------|---------------------------|------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661040099873 | Warranty (months) | 12 |
| Country | Germany | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 1.16 |
| Bare Unit / Weight (kg) | 1 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Light blue | Shoe Size | 38 |

 28/04/2025 22:54:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr