

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

medi Plogs

Ref. WSAM_39V



📅 28/04/2025 22:39:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31
✉ info@gsh-med.fr

MEDILOGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_39V**

=====


Description


- Clogs suitable for modern hospitals, operating theatres and intensive care units.
- Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
 - The back of the heel is sloped and rounded to make walking easier
 - Ventilation system allows air to circulate
 - Toe prints to enhance hold and spurs to aid the circulation
 - 100% CFC-free polyurethane
 - Washable by hand or machine at 60°C
 - Sizes: 35 - 47

Features

CE Marking	NONE	Medical device	No
Bar code	3661040099972	HS Code	64029993
Warranty (months)	12	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.16
Bare Unit / Weight (kg)	1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Green	Shoe Size	39

 28/04/2025 22:39:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr