

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

medi Plogs

Ref. WSAM\_42V



📅 28/04/2025 22:36:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)

# MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM\_42V**

=====


## Description


- Clogs suitable for modern hospitals, operating theatres and intensive care units.
- Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
  - The back of the heel is sloped and rounded to make walking easier
  - Ventilation system allows air to circulate
  - Toe prints to enhance hold and spurs to aid the circulation
  - 100% CFC-free polyurethane
  - Washable by hand or machine at 60°C
  - Sizes: 35 - 47

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040100159	HS Code	64029993
Warranty (months)	12	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.16
Bare Unit / Weight (kg)	1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Green	Shoe Size	42


28/04/2025 22:36:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)