

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

mediPlogs

Ref. WSAM_43BC



28/04/2025 22:46:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_43BC**

=====

Description

- Clogs suitable for modern hospitals, operating theatres and intensive care units.
- Products designed and manufactured in collaboration with researchers and doctors.
- Antistatic clogs with non-slip and non-chafing soles
 - The back of the heel is sloped and rounded to make walking easier
 - Ventilation system allows air to circulate
 - Toe prints to enhance hold and spurs to aid the circulation
 - 100% CFC-free polyurethane
 - Washable by hand or machine at 60°C
 - Sizes: 35 - 47

Features

CE Marking	NONE	Medical device	No
Bar code	3661040100173	Warranty (months)	12
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.16
Bare Unit / Weight (kg)	1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Light blue	Shoe Size	43

 28/04/2025 22:46:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr