

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

mediPlogs

Ref. WSAM_44BC



28/04/2025 22:10:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE



Ref. **WSAM_44BC**

=====

Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

Features

| | | | |
|------------------------------------|---------------|---------------------------|------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661040100234 | Warranty (months) | 12 |
| Country | Germany | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 1.16 |
| Bare Unit / Weight (kg) | 1 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Light blue | Shoe Size | 44 |

28/04/2025 22:10:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr