


CARDIOTOCOGRAPHE C11 WITH TWIN OPTION AND VCT

COMEN

Ref. BCTOCO12



CE
1639



DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
3 ANS
3
GUARANTEE



03/05/2025 02:16:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CARDIOTOCOGRAPHE C11 WITH TWIN OPTION AND VCT

COMEN

Ref. BCTOC012

A cardiotocograph very easily transportable with its compact design and elegant mini Eiffel carrying bag!

Description

The C11 is a fetal and maternal monitor that performs fetal monitoring monitoring. Thus it measures the cadence of the heart of the fetus as well as the contractions of the uterus of the mother.

It stands out for its portability: compact, lightweight and easily transportable with its handle or elegant Eiffel mini bag (optional). It is perfectly suited to private midwives, gynaecology/obstetrics practices, as well as maternity services, thanks to its trolley with wheels that allows it to be easily moved to the bedside of the future mother.

Compact design with an elegant optional mini Eiffel bag

5.7" color LCD screen adjustable from 0 to 90°

Automatic recognition of different probes to avoid malfunction

Very sensitive (9 crystals) and waterproof (IP67) probes, very useful at the time of delivery

Fetal motion detector (automatic and manual)

Fast and adjustable thermal printer

VCT calculation and advanced NST/Fischer/Fischer scores/ Krebs

Dimension: 28.7 x 23.2 x 9 cm Weight: 2.2kg

Warranty: 3 years

Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219233	HS Code	9018193010
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	29	Bare Unit / Width (cm)	23.2
Packed unit / Height (cm)	38	Bare Unit / Height (cm)	9
Packed unit / Length (cm)	45	Bare Unit / Length (cm)	28.7
Packed unit / Weight	4.55	Bare Unit / Weight (kg)	2.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 03/05/2025 02:16:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr