

ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS

Ref. L35587



06/05/2025 23:57:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS



Ref. L35587

Description

Adult tips for SPENGLER peak flow meter.
Dimensions of the tips 30mm (diameter) by 67 mm (length).
Sold in boxes of 100 pieces.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446001681	HS Code	48229000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.83
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/05/2025 23:57:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr