

## MINI BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

Ref. MDOUES02

*De Boissy*  
— DEPUIS 1982



07/05/2025 00:03:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MINI BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

*De Boissy*  
— DEPUIS 1982

Ref. MDOUES02

The mini version of the OUessant Medical Bag

### Description

Eco-responsible, customisable and 'super-organised' for a trendy look with all the practical, modern options.

All the fabrics, both exterior and interior, are Global Recycled Standard (GRS) certified, made from recycled materials and remain washable, hard-wearing and water-repellent.

Its deep blue, inspired by the shores of Ushant, is enhanced by black grained leather-effect finishes, giving it an elegant, professional look.

While De Boissy was keen to pay tribute to carers and take inspiration for its bags from the shape of the traditional doctors' leather satchel, the look of the Ouessant bags contrasts sharply with the usual classicism. De Boissy has opted for elegance, trendiness and modernity. And the fact remains that they have been designed to be 'ultra practical', to be the best ally of healthcare professionals during their medical visits at home or in the community.

- The mini version of the Etretat and Eiffel Medical Bags
- The main space inside also includes 6 pockets, a divider and 2 mobile compartments: 1 transparent and 1 isothermal, supplied with its reusable cold gel pouch.
- The belly pocket measures 31x15 cm
- Still on the outside: the 3 small zipped pockets, plus the one with elastic for the gel bottle, the detachable-customisable grained leather-look pouch, and the external USB port with integrated charging cable
- Same handles, shoulder strap, name tag and non-slip bottom as the large format
- Dimensions (W x H x D): 32 x 28 x 20 cm
- Weight: 1.6 kg
- Maximum recommended contents: 4 kg / 15 litres
- Sold empty

The De Boissy bag comes in a dust bag to protect it when stored. Grey with patterns, this accessory is made from recycled materials with a drawstring fastening.

📅 07/05/2025 00:03:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)


☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)

**MINI BAG OUessant, DE BOISSY,  
WITH BLACK LEATHER POCKET***De Boissy*  
— DEPUIS 1982Ref. **MDOUES02****Features**

CE Marking	NONE	Medical device	No
Bar code	3700446042424	HS Code	4202220000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	54	Packed unit / Height (cm)	42
Packed unit / Length (cm)	57	Packed unit / Weight	1.05
Bare Unit / Weight (kg)	1.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

 07/05/2025 00:03:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)