


CONVERTIBLE BACK PACK OUESSANT DE BOISSY, WITH BLACK LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. **MDOUES03**



06/05/2025 23:49:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CONVERTIBLE BACK PACK OUESSANT DE BOISSY, WITH BLACK LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDOUES03

It can be carried as a rucksack, in the hand or on the shoulder - a real Swiss Army knife!

Description

Eco-responsible, customisable and 'super-organised' for a trendy look with all the practical, modern options.

All the fabrics, both exterior and interior, are Global Recycled Standard (GRS) certified, made from recycled materials and remain washable, hard-wearing and water-repellent.

Its deep blue colour, inspired by the shores of Ushant, is enhanced by black grained leather-effect finishes for an elegant, professional look.

OUESSANT bags have been designed to be 'ultra practical', to be the best ally of healthcare professionals during their medical visits at home or in the community.

- Its main compartment is modular and contains 4 zipped, removable compartments (secured with Velcro): 3 transparent for quick identification of contents, and 1 insulated, supplied with its reusable cold gel pouch, for up to 20 ampoules.
- An extension zip can be used to enlarge the main compartment to hold more equipment
- A second padded area for your computer/tablet and documents
- This model also has a large belly pocket (45x21 cm) that opens completely, with 2 zipped storage nets and 2 flat pockets on one side, and several elastic loops on the other side for storing small medical supplies: tweezers, scissors, penlight, thermometer, bandages, pens, etc...
- Two external zipped pockets, plus one with an elastic band for a bottle of hydro-alcoholic gel, a detachable and customisable grained leather-look pouch (fancy printed pouch sold separately)
- An external USB port with integrated charging cable means you can plug in a back-up battery and never run out of phone or forget your charger at a patient's bedside.
- Padded shoulder straps, handles and shoulder strap are sturdy and concealable, so you can wear it in the way that suits you best
- Dimensions (W x H x D): 30 x 45 x 15/19 cm
- Weight: 2.2 kg
- Maximum recommended content: 6 kg / 19-25 litres

The De Boissy bag comes in a dust bag to protect it when it's put away. Grey with patterns, this accessory is made from recycled materials and has a drawstring fastening.

📅 06/05/2025 23:49:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31
✉ info@gsh-med.fr

CONVERTIBLE BACK PACK OUESSANT DE BOISSY, WITH BLACK LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDOUES03

Features

| | | | |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446042431 | HS Code | 4202129000 |
| Warranty (months) | 12 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 50 | Packed unit / Height (cm) | 35 |
| Packed unit / Length (cm) | 50 | Packed unit / Weight | 1.86 |
| Bare Unit / Weight (kg) | 1.86 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |

06/05/2025 23:49:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr