


## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519011



07/05/2025 21:59:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519011

=====

### Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

#### CUFF

- Cotton CUFF, velcro fastening
- Size XXS (10-15 cm), XS (15-22 cm), S (20-28 cm, ref : 519029), M (26-34 cm, ref 519012), or L (32-42 cm), S-M-L (ref :519050)
- Available in grey, navy, red or patterns
- Supplied with carrying case

### Features

|                   |               |                |                               |
|-------------------|---------------|----------------|-------------------------------|
| CE Marking        | CE0459        | Medical device | Class Im (measuring function) |
| Bar code          | 3700446007218 | HS Code        | 90189010                      |
| Warranty (months) | 36            | Country        | France                        |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.706 |

#### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 70 |
| Minimum humidity    | 85  | Maximum humidity | 85 |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

07/05/2025 21:59:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)