

MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX



Ref. WMBOU07_94



13/05/2025 23:34:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX



Ref. **WMBOU07_94**

The eco-responsible men's medical tunic

Description

- White men's tunic with short sleeves,
- Eco-responsible: 50% Lyocell (eucalyptus fiber) and 50% Polyester,
- Officer collar,
- Twill, 215 g / m2,
- 3 plated pockets,
- Buttoning by pressure.
- Available in :
Size 0 (WMBOU07_90)
Size 1 (WMBOU07_91)
Size 2 (WMBOU07_92)
Size 3 (WMBOU07_93)
Size 4 (WMBOU07_94)
Size 5 (WMBOU07_95)

The Holtex eco-friendly and responsible medical tunic brings extra comfort and strength to a classic polycotton tunic: superior quality, increased strength, easy maintenance and maximum ease. Its cut is traditional and modern at the same time. Its officer collar, its short sleeves mounted lapels and its 3 plated pockets are nicely finished with a white border which makes this men's tunic a very neat garment.

Lyocell is a natural material, made from eucalyptus wood pulp.

- Soft, lightweight, comfortable and slightly elastic,
- Easy to maintain, wrinkle resistant, it does not lint.
- Antibacterial, hypoallergenic, breathable, it absorbs moisture (up to 50% more than cotton): it is ideal for sensitive skin.
- Environmentally friendly: Eucalyptus cultivation requires much less water and pesticides than cotton.

Features

| | | | |
|------------|---------------|----------------|------------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661041009512 | HS Code | 6211433100 |
| Country | Tunisia | | |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | | | |
|------------------------|----------------------|----------|--------------|
| Bleaching | Bleaching prohibited | Washing | Wash at 60°C |
| Model | Male model | Cleaning | Dry cleaning |
| Packed unit / Quantity | 1 | Size | Size 4 |

10/03/2023 23:34:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr