

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

CE  
0459



27/05/2025 19:31:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

=====

### Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

### Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446007287 | HS Code                   | 90189010                      |
| Warranty (months)           | 60            | Country                   | France                        |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | 14            | Packed unit / Height (cm) | 8                             |
| Packed unit / Length (cm)   | 20            | Packed unit / Weight      | 0.75                          |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |
| Minimum humidity            | 20            | Maximum humidity          | 85                            |
| OTHER FEATURES              |               |                           |                               |
| Packed unit / Quantity      | 1             |                           |                               |

27/05/2025 19:31:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)