

Ref. E1.417.012



26/05/2025 19:18:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

=====

Ref. E1.417.012

LUXAMED® laryngoscope, the quality for all applications!

Description

Laryngoscopes are used for intubation in operating room and in every emergency cases such as respiratory emergency and resuscitation. They are also used to perform laryngoscopies for therapeutic purpose.

The main part of the LUXAMED® laryngoscope is the handle. Its quality is especially characterized by illumination and longevity.


Lateral grooves ensure easy cleaning and secure grip in the hand.

They are available with different types of blades: Macintosh, Miller, and McCoy with different sizes.

Features

| | | | |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 4260381960348 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -10 | Max temperature | 45 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |

📅 26/05/2025 19:18:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr