

## DDM ELECTRONIC PNEUMATIC TOURNIQUET EASY PUMP® DUAL 2 PRESSURE CIRCUITS

Ref. **G10904**

**DDM**

CE  
0459

FABRICATION  
FRANÇAISE

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**2** ANS  
2 YEARS  
GUARANTEE



27/05/2025 19:41:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## DDM ELECTRONIC PNEUMATIC TOURNIQUET EASY PUMP® DUAL 2 PRESSURE CIRCUITS



Ref. G10904

The only electronic tourniquet manufactured in France!

### Description

The EASY PUMP® DUAL tourniquet is DDM's new generation of electronic tourniquet providing a safe, bloodless operating field for lower limb surgery. A model with 2 independent pneumatic circuits to be used with a single or double cuff or with 2 single cuffs. It allows 2 cuffs to be inflated at different pressures for bilateral surgery or IntraVeinuous Regional Anesthesia (IVRA).

Easy to use and intuitive:

- Pre-settings of time, pressure, languages...
- Large screen 7" with IPS touch screen technology (1024 x 600 pixels) with no viewing angle
- Intense brightness (700cd/m²)
- Highly intuitive icons for setup and use.
- Distinct color codes for different circuits

Safe for the patient:

- Integrated LOP function without SpO2 sensor, with corrective adjustment.
- Autotest upon switch on
- Time and pressure preset by default and can be modified during operation.
- IVRA mode and possible guidance function.
- Leakage compensation.
- Visual and audible alarms.

### Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700468214960	HS Code	9018902090
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	19
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	30
Packed unit / Weight	nc	Bare Unit / Weight (kg)	5.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	20	Max temperature	60

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)