



MASTER PALM 2 - HANDHELD OXIMETER WITH INTERCHANGEABLE SENSORS



Ref. 223000



01/06/2025 21:04:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MASTER PALM 2 - HANDHELD OXIMETER WITH INTERCHANGEABLE SENSORS



Ref. 223000

The MASTER PALM 2® pulse oximeter is the pulse oximeter specially designed for professionals

Description

The pulse oximeter Master Palm 2® is suitable for health professionals use.

Besides its lightness and exceptional battery life, you will enjoy the speed and accuracy of its readings and its versatility.

Its technology allows a quick intervention on all types of patients from adult to neonatal, to precisely check their blood oxygen saturation and pulse rate.

- High definition screen, with auto-rotating function (activated by gravity), intend to optimize ergonomics by providing a vertical or horizontal layout.

- 2 display modes: Numbers and SpO2 curve. Display pulse, plethysmography curve, SpO2 and SpO2 curve.

- Fully configurable: Set the alarm, display, sound, system settings.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446025380	HS Code	9018199000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.718
Bare Unit / Weight (kg)	0.718		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

01/06/2025 21:04:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr