

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/LARGE ADULT (XL)

Ref. 337870



01/06/2025 21:07:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/LARGE ADULT (XL)



Ref. 337870

Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Welch Allyn, BCI (minition plus) patient monitors and Suntech 247®

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028299	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.134
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

01/06/2025 21:07:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr