

TARNIER OBSTRETRICAL FORCEPS BOX



Ref. IBO109



03/06/2025 00:26:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

TARNIER OBSTRETRICAL FORCEPS BOX



Ref. IBO109

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040081663	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	20.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	41
Packed unit / Weight	2.849		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:26:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr