

SPO2 SENSOR PAEDIATRIC CLIP FOR COMEN NC3/NC5/C31

COMEN

Ref. DBCPEDI01



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
1 YEAR
1
GUARANTEE

02/06/2025 21:51:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPO2 SENSOR PAEDIATRIC CLIP FOR COMEN NC3/NC5/C31

COMEN

Ref. DBCPEDI01

Designed for accurate and fast measurement of oxygen saturation

Description

Technical specifications: % oxygen 90-100% accuracy +/- 2% 70-90%

accuracy +/- 3%
<70% not specified

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219424	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 21:51:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr