

ADULT SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C32

COMEN

Ref. DBSADUL01



CE
0537

DISPOSITIF MEDICAL
CLASSE
II^b
MEDICAL DEVICE

GARANTIE
1
1 YEAR
GUARANTEE

02/06/2025 23:50:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ADULT SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C32

COMEN

Ref. DBSADUL01

Designed to ensure accurate and fast measurement of oxygen saturation.

Description

Technical specifications: % oxygen: 90-100% Accuracy +/-2% 70-

90% Accuracy +/-3%
<70% not specified.

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219431	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 23:50:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr