


PEDIATRIC SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C33

COMEN

Ref. DBSPEDI01



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE

03/06/2025 00:30:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PEDIATRIC SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C33

COMEN

Ref. DBSPEDI01

Designed to ensure accurate and fast measurement of oxygen saturation.

Description

Technical specifications% oxygen 90-100% Accuracy +/- 2% 70-90%

Accuracy +/-3%
<70% not specified.

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219448	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 03/06/2025 00:30:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr