

TEMPO II ADULT CUFF (26 TO 35 CM)

Ref. DS99645



03/06/2025 00:31:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

TEMPO II ADULT CUFF (26 TO 35 CM)



Ref. DS99645

Immersible and decontaminable armband

Description

Velcro Cuff, without pocket, submersible, with cap

Adult cuff (26 to 35 cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040062242	HS Code	9018901000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.093
Bare Unit / Weight (kg)	0.093		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:31:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr