


PERSONAL SCALES WITH COLUMN FOR CONVENIENT, ERGONOMICALLY-OPTIMISED (MAX. 250KG)

KERN®

Ref. EKMPK01



03/06/2025 00:52:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

PERSONAL SCALES WITH COLUMN FOR CONVENIENT, ERGONOMICALLY-OPTIMISED (MAX. 250KG)

KERN®

Ref. EKMPK01

Ideal measuring instrument for elderly care, paramedics, fitness clubs or in the private sector and can handle weighing overweight persons weighing up to 250 kg

Description

- The display device is easy to operate through the column which is firmly fixed to the scale and the weight can be read easily, even under poor lighting conditions
- Sturdy glass weighing plate for secure, comfortable standing position ; Easy and hygienic cleaning
- Large high-contrast display that is easy to read (height of digits 28 mm)
- Particularly flat design facilitates access
- Easy to operate: simply tap on it – stand on it – read the weight
- Integrated AUTO-OFF function to preserve the batteries
- Total dimensions (including column) L × W × H 350 × 420 × 955 mm
- Batteries included, 4 × 1.5V AA, service life up to 20 h
- Net weight approx. 3.8 kg

Features

CE Marking	CE	Medical device	Class Im (measuring function)
Bar code	3661040195414	HS Code	84231090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	42
Packed unit / Height (cm)	45	Bare Unit / Height (cm)	95.5
Packed unit / Length (cm)	102	Bare Unit / Length (cm)	35
Packed unit / Weight	6.4	Bare Unit / Weight (kg)	3.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:52:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr