

BINKHORST CATARACT ASPIRATION CANNULA STRAIGHT 0.5 MM

Ref. IC92106



02/06/2025 22:56:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

BINKHORST CATARACT ASPIRATION CANNULA STRAIGHT 0.5 MM



Ref. IC92106

Features

CE Marking	CE	Medical device	Class IIa
Bar code	3661040013114	HS Code	9018.9090
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	5.7
Packed unit / Weight	0.033	Bare Unit / Weight (kg)	0.033
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 22:56:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr