

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M10500

DDM

CE



03/06/2025 23:12:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM



Ref. M10500

The single use pressure cuff made in France !

Description

The pressure infusion cuffs are medical devices designed for daily use in operating theaters, Intensive Care units or emergency transport vehicles. They can be used to manually accelerate infusion of blood and solutions when gravity alone cannot achieve the desired flow rates.

Easyfuse pressure cuffs are single use and compatible with different sizes of blood/solution bags: 500cc, 1000cc and 3000cc.

Optimized monitoring

- **Single use**, no cleaning and no disinfection
- blood/solution bag clearly visible
- pressure indicator with automatic overpressure regulator
- Controlled decompression facilitated by 3 ways valve
- Resistant as made of TPU

All our pressure infusion cuffs are made in France in our Castelsulier (Agen) manufacturing site..

 03/06/2025 23:12:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM





Ref. M10500

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203810	HS Code	98109084
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			


03/06/2025 23:12:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr