


SPO2 SENSOR PEDIATRIC FINGER CLIP FOR MASTER PALM 2



Ref. 222251



03/06/2025 00:48:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPO2 SENSOR PEDIATRIC FINGER CLIP FOR MASTER PALM 2



Ref. 222251

Description

Rigid pediatric SpO2 sensor for MASTER PALM 2 oximeter

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030216	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.093
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


03/06/2025 00:48:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr