

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC SILICON FOR MP3




Ref. 222265



CE
0537

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
1
1 YEAR
GUARANTEE



02/06/2025 23:30:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC SILICON FOR MP3



Ref. 222265

Description

NELCORR OXIMAX compatible pediatric silicon spo2 sensor for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030186	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.084
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 23:30:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr