

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520621



CE
0459

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
3
3 YEARS
GUARANTEE



03/06/2025 00:06:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520621

=====

Description

HIGN ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: XS (10-15 cm), XS (15-22 cm), S (20-28 cm), M (26-34 cm), L (32-42 cm) ou XL (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008406	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.667

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
---------------------	-----	-----------------	----

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr