

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520630



CE  
0459



02/06/2025 23:16:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520630

=====

### Description

#### HIGN ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

#### SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

#### COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

#### CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: XS (10-15 cm), XS (15-22 cm), S (20-28 cm), M (26-34 cm), L (32-42 cm) ou XL (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

### Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446007249 | HS Code                   | 90189010                      |
| Warranty (months)           | 36            | Country                   | France                        |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 1.122                         |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)