

LIAN METAL / GREY NYLON CUFF / OBESE (XL)



Ref. 521312



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3
3 YEARS
GUARANTEE

03/06/2025 01:09:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN METAL / GREY NYLON CUFF / OBESE (XL)



Ref. 521312

=====

Description

Large dial, ambidextrous model (simply rotate to change hands).Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

Available in :

- Adult (M), cotton, velcro, red (521317)
- Adult (M), cotton, velcro, grey (521311)
- Adult large (L), cotton, velcro, grey (521310)
- Child/Small Adult (S), nylon, velcro, grey (521316)
- Multibrassard (S, M, L) velcro (521330)

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446006631	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.551
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 03/06/2025 01:09:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr