

BLACK NYLON SOFT BAG SPENGLER, MEDIUM

Ref. 535020



03/06/2025 00:53:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BLACK NYLON SOFT BAG SPENGLER, MEDIUM



Ref. 535020

Description

Nylon pouch 18x13x6 cm with eclair® zip closure

Features

CE Marking	NONE	Medical device	No
Bar code	3700446032395	HS Code	90330000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:53:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr