

LARGE ADULT CUFF 27-42CM FOR SPOT-CHECK 2VS



Ref. 546020



03/06/2025 00:00:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LARGE ADULT CUFF 27-42CM FOR SPOT-CHECK 2VS



Ref. 546020

Description

Adult wide armbands (27-42cm) for spot-check station 2VS from SPENGLER.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029678	HS Code	9018901000
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.091
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


03/06/2025 00:00:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr