

PEDIATRIC CUFF 18-26 CM FOR 4VS MONITOR

Ref. 551621



02/06/2025 23:47:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PEDIATRIC CUFF 18-26 CM FOR 4VS MONITOR



Ref. 551621

Description

Pediatric cuff for PNI module of SPENGLER 4VS vital signs monitor.
Dimensions (18-26cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446033057	HS Code	9018901000
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.049
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 23:47:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr