

SPO2 ADULT SENSOR FOR 4VS MONITOR



Ref. 551630



CE
0537



02/06/2025 22:08:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPO2 ADULT SENSOR FOR 4VS MONITOR



Ref. 551630

Description

SPO2Adult Rigid Sensor for SPENGLER 4 VS Vital Signs Monitor

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031176	HS Code	90181990
Warranty (months)	6	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.103
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		


02/06/2025 22:08:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr