

SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR



Ref. 551635



CE
0537



03/06/2025 00:21:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPO2 PEDIATRIC SILICON SENSOR FOR 4VS MONITOR



Ref. 551635

Description

SPENGLER Pediatric Silicone SPO2 Sensor for Vital Signs Monitor 4 VS

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446031220	HS Code	90181990
Country	Uruguay		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.083
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-40	Max temperature	55
Minimum humidity	5	Maximum humidity	95
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:21:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr