

## EXTENSIN CABLE FOR IOP SENSOR (2M)



Ref. 554111



03/06/2025 00:43:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## EXTENSIN CABLE FOR IOP SENSOR (2M)



Ref. 554111

### Features

|                             |                |                           |              |
|-----------------------------|----------------|---------------------------|--------------|
| CE Marking                  | NONE           | Medical device            | No           |
| Bar code                    | 28809315676540 | HS Code                   | 9018 90 2090 |
| Warranty (months)           | 12             |                           |              |
| DIMENSIONS AND WEIGHT       |                |                           |              |
| Packed unit / Width (cm)    | nc             | Packed unit / Height (cm) | nc           |
| Packed unit / Length (cm)   | nc             | Bare Unit / Length (cm)   | 2            |
| Packed unit / Weight        | 0.2            |                           |              |
| SPECIFIC STORAGE CONDITIONS |                |                           |              |
| Minimum temperature         | NA             | Max temperature           | NA           |
| Minimum humidity            | NA             | Maximum humidity          | NA           |
| OTHER FEATURES              |                |                           |              |
| Packed unit / Quantity      | 1              |                           |              |

03/06/2025 00:43:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)