


## SMARTLIGHT OTOSCOPE BLACK

Ref. 570315



02/06/2025 21:51:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SMARTLIGHT OTOSCOPE BLACK



Ref. 570315

### The otoscope with conventional lighting

### Description

Conventional lighting. Lightweight and robust, this otoscope is perfect for ENT auscultations especially on the move. The XHL lamp diffuses a concentrated and homogenous light beam (optional LED lamp available). Fit with an ABS and Brass head, SMARTLIGHT also comes with end clamp and built-in switch, LR6 AA batteries included.

### Features

CE Marking	CE	Medical device	No
Bar code	3700446041557	HS Code	90189042
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Packed unit / Height (cm)	3
Packed unit / Length (cm)	10	Packed unit / Weight	0.23
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 21:51:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)