

## LITESTICK PENLIGHT SPENGLER

Ref. 570622



02/06/2025 22:03:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LITESTICK PENLIGHT SPENGLER



Ref. 570622

### Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

### Features

|                             |               |                           |        |
|-----------------------------|---------------|---------------------------|--------|
| CE Marking                  | CE            | Medical device            | No     |
| Bar code                    | 3700446009878 | HS Code                   | 901890 |
| Warranty (months)           | 24            |                           |        |
| DIMENSIONS AND WEIGHT       |               |                           |        |
| Packed unit / Width (cm)    | 5             | Packed unit / Height (cm) | 5      |
| Packed unit / Length (cm)   | 20            | Packed unit / Weight      | 0.046  |
| SPECIFIC STORAGE CONDITIONS |               |                           |        |
| Minimum temperature         | 10            | Max temperature           | 40     |
| OTHER FEATURES              |               |                           |        |
| Colors                      | Green         | Packed unit / Quantity    | 1      |

02/06/2025 22:03:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)