



IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_N



03/06/2025 01:03:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

IDEAL+ DOUBLE HEAD ADULT, BLUE



Ref. DH7073_N

=====

Description

- Stethoscope with double metal horn adult (O: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073_B), Black (DH7073_N), Red (DH7073_R) or Pink (DH7073_RC) for the paediatric version (O: 35/30 mm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000083	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.35
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

03/06/2025 01:03:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr