

SPENGLER ARMBAND, VAQUEZ- LAUBRY NANO POCKET



Ref. 330220

 02/06/2025 22:44:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPENGLER ARMBAND, VAQUEZ-LAUBRY NANO POCKET



Ref. 330220

Description

Made of black Oxford-grade nylon, with inner bladder and dual tubing (10 cm & 60 cm), with split leather manometer holder, and Velcro fastening.

Child/Small adult (S), nylon, velcro, black : 330220

Adult (M), nylon, velcro, black : 330210

Adult (L), nylon, velcro, black : 330200

Features

CE Marking	NONE	Medical device	No
Bar code	3700446028565	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 22:44:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr