

**LIAN CUFF/1
TUBE/NYLON/BLACK/LARGE ADULT
(L)**



Ref. 333000



 03/06/2025 00:39:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIAN CUFF/1 TUBE/NYLON/BLACK/LARGE ADULT (L)



Ref. 333000

Description

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color.

Washable at 30-40° C.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446008765	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.21
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/06/2025 00:39:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr