

HEATING BLANKET


Ref. CB48-3+BCB48

CE
0459

FABRICATION
FRANÇAISE

DISPOSITIF MEDICAL
CLASSE
I Ib
MEDICAL DEVICE

GARANTIE
1 YEAR
GUARANTEE



03/06/2025 23:11:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

HEATING BLANKET

Ref. CB48-3+BCB48

The patient warming system for perioperative care

Description

Because studies have shown that intraoperative hypothermia can lead to complications for patients. GSH offers you ECM brand warming blankets, which are active patient warming systems perfectly suited to the various perioperative phases.

Eco-responsible

Adapted technology: circuit integrated into a cartridge, no pulsed air, thus preventing the presence of particles in the operating environment.

Inertia when the power is cut off, keeping the patient warm for a few more minutes.

Reusable: reduces consumable waste and is more economical.

Ease of use

Reusable

Several sizes to suit every patient

Select the temperature from 30 to 40° using the control knob on the unit and switch on the unit.

Covered with a PU cover for easy cleaning

In complete safety

Fitted with 8 independent temperature sensors for precise operation and maximum safety

with visual and audible alarms

Low-voltage power supply to prevent interference with monitors

Made in France

We have been designing and manufacturing patient warming systems for hospitals since 1991, at our factory in Saint Leonard de Noblat (Limoges).

Available in 3 sizes:

-half-body/paediatric (Ref : CB48-2+BCB48)

-upper body (Ref: CB48-3+BCB48)

-full body (ref :CB48-11+BCB48)

Translated with DeepL.com (free version)

📅 03/06/2025 23:11:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

HEATING BLANKET

Ref. CB48-3+BCB48

Features

CE Marking	CE0459	Medical device	Class IIb
HS Code	90189084	Warranty (months)	12
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	50	Packed unit / Height (cm)	11
Packed unit / Length (cm)	150	Packed unit / Weight	nc
Bare Unit / Weight (kg)	3.4		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	10	Max temperature	40
Minimum humidity	30	Maximum humidity	90
OTHER FEATURES			

03/06/2025 23:11:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr