

CHARGING BASE STATION FOR THERMOSCAN PRO 6000 TYMPANIC THERMOMETER

WelchAllyn®

Advancing Frontline Care™

Ref. DWPRO6K07



02/06/2025 22:14:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CHARGING BASE STATION FOR THERMOSCAN PRO 6000 TYMPANIC THERMOMETER

WelchAllyn®

Advancing Frontline Care™

Ref. DWPRO6K07

Switch from Alkaline Batteries Mode to Rechargeable Battery

Description

- Allows you to automatically recharge the thermometer thanks to a rechargeable battery.
- Offers storage space for 200 probe protectors
- Has an electronic security function – to configure: if the thermometer is not returned to the charging station within the selected allotted time, the safety function locks it
- Can be fixed in wall bracket
- Comes with a rechargeable battery

Features

CE Marking	NONE	Medical device	No
Bar code	00732094325324	HS Code	850410
Warranty (months)	36	Country	Mexico
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	9
Packed unit / Length (cm)	25	Packed unit / Weight	0.7
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 22:14:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr