

AUTOCLAVABLE DOUBLE CUFF ARM SILICON BLADDER ADULT 67 X 17 CM

Ref. GHP552812

HpmTM

CE



 03/06/2025 00:04:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AUTOCLAVABLE DOUBLE CUFF ARM SILICON BLADDER ADULT 67 X 17 CM



Ref. GHP552812

Operating room equipment: double pocket autoclavable cuff for two-way electronic pneumatic tourniquet. Sold individually.

Description

Nylon cuff with double silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilization.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446033767	HS Code	9018 90 2090
Warranty (months)	12	Country	Sweden
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.51
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


 03/06/2025 00:04:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr