

SAUTTER HYDRODISSECTION CANNULA

Ref. IC92118



02/06/2025 23:49:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SAUTTER HYDRODISSECTION CANNULA



Ref. IC92118

Features

CE Marking	CE	Medical device	Class IIa
Bar code	3661040145006	HS Code	901850
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Bare Unit / Width (cm)	6.7
Packed unit / Height (cm)	1	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20.2
Packed unit / Weight	0.014	Bare Unit / Weight (kg)	0.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 23:49:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr