

## ELITE BAGS-EXTREME- POLYESTER/WATERPROOF-RED



Ref. MEEXTR02\_R



02/06/2025 22:17:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ELITE BAGS-EXTREME-POLYESTER/WATERPROOF-RED



Ref. MEEXTR02\_R

### Description

For users of emergency services specialized in first aid.

Interior totally adaptable to the needs of the user thanks to extractable velcro separators.

Elastic bands in all interior flat surfaces.

Reflective strips night and day. 2 external pockets of medium size, 2 external pockets of large capacity, with elastic bands.

4 extractable compartments with transparent frontal part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: various surgical equipment, yellow: pediatric equipment).

Isothermal amp maker with 27 elastics.

Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana.

Material: 1000D polyamide fabric

Colour: Red.

Weight: 2.86 kg

Outdoor measures: 47 x 22 x 28 cm / 23.23 l. Maximum recommended capacity: 12 Kg

### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540300	Warranty (months)	24
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Bare Unit / Width (cm)	22
Packed unit / Height (cm)	39	Bare Unit / Height (cm)	28
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	47
Packed unit / Weight	3.8	Bare Unit / Weight (kg)	2.86
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/06/2025 22:17:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)