

CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)




Ref. 260051

CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
5
5 YEARS
GUARANTEE



06/06/2025 21:42:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CARDIOMATE EVI 10" SPENGLER ECG TABLET (WITH PRINTER BASE)



Ref. 260051

Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 10" tablet, a product of a new technology, totally autonomous.

Description

A very great simplicity of use: - 10 "color touch screen- Portable: compact size and light weight, easy to carry: 1.2 kg and 272.6 x 181.2 x 25.8 mm)- Large capacity battery, rechargeable on mains or on the recording base- Allows to print ECGs directly on a printer via wifi network or on the recorder base High-performance technology:

High accuracy: sampling speed 32,000 Hz

- Detection and marking of the

pacemaker

- Automatic interpretation (Glasgow) and calculation of intervals- Input of

patient information using a keyboard or barcode reader-

Visualization of the 12 tracks (15/18 other tablet optional) at the same time as automatic interpretation and interval calculations.

- Data export via wifi or USB port under .pdf or .xml

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6971049341796	HS Code	9018110000
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	24	Packed unit / Height (cm)	36
Packed unit / Length (cm)	36	Packed unit / Weight	7
Bare Unit / Weight (kg)	3.3		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES

06/06/2025 21:42:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr